



# Hitch

## Compressible Bone Matrix

Fill bone voids with the Hitch Compressible Bone Matrix.



### KEY BENEFITS

- 100% demineralized cancellous bone
- Can compress to half its size, then re-expand to fill bone voids
- Strips, cubes and chips available for a variety of applications
- Decreases operating room time
- Excellent osteoconductivity
- Supports rapid vascularization and new bone growth

The Hitch Compressible Bone Matrix is an expansive material composed of 100% cancellous bone. Its ability to compress to half its size, then re-expand to fill tight-space bone voids, makes it a versatile tool for use throughout the skeletal system. Upon re-expanding after placement, the sponge allows for better fit and reduces potential for migration. The conductivity supports vascularization and stimulation of new bone growth.

STYLES	SIZES	REORDER NUMBER
Cubes	10mm	TTH101010
	14mm	TTH141414
Strips	20x50x5mm	TTH205005
Chips	4-10mm 10cc	TTH010410
	4-10mm 15cc	TTH015410

Learn more about how Hitch significantly improves the healing process—and patient outcomes—by calling **888.494.4441** or visiting **tidesmedical.com**.



tidesmedical

888.494.4441 [tidesmedical.com](http://tidesmedical.com)